

## Pharmacological Characterization Utilizing Manual Patch Clamp

Independent analysis of the **spontaneous and pacing stability** of vCor.4U™ and their **predictive response to selective pharmacological agents**



## OVERVIEW

Extensive cardiac ion channel profiling at Metrion using selective pharmacological tools and manual patch clamp has enabled us to define the electrophysiological characteristics of vCor.4U™ human iPSC-derived cardiomyocytes (iPSC-CM).

vCor.4U™ elicit spontaneous action potentials (AP) and can be paced at a range of frequencies (0.2 - 1 Hz; Figure 1). The electrophysiological characteristics of evoked AP resemble primary human cardiomyocytes, with the cells having an average resting membrane potential of -73 mV and action potential duration at 90 % of repolarization (APD90) of 426 ms when paced at 1 Hz (Figure 1C).

vCor.4U™ AP are stable over time in both spontaneous and evoked recording conditions confirming the suitability of vCor.4U™ for compound screening experiments (Figure 2). Additionally, we confirm that vCor.4U™ express a predominantly ventricular phenotype and exhibit the appropriate pharmacology in response to core cardiac channel modulators, including early after depolarizations (EADs) following  $I_{Kr}$  inhibition (Figure 3).

In the light of the incoming CiPA guidelines, Axiogenesis vCor.4U™ iPSC-CM represent a suitable model for *in vitro* pre-clinical cardiac safety evaluation of compounds to identify potential human proarrhythmic liabilities.

## vCOR.4U™ HAVE ELECTROPHYSIOLOGICAL CHARACTERISTICS THAT RESEMBLE PRIMARY HUMAN CARDIOMYOCYTES




Figure 1: Characteristics of spontaneous and evoked vCor.4U™ action potentials

Representative traces of spontaneous (A) and evoked (B) AP (1 Hz) recorded from vCor.4U™ iPSC-CM under control conditions. (C) Average parameters for spontaneous (n = 89) and evoked (n = 75) AP in control conditions. RMP, resting membrane potential;  $V_{max}$ , maximum upstroke velocity; APA, action potential amplitude; APD20-90, action potential duration at 20, 50 or 90 % repolarization.

## vCOR.4U™ ADVANTAGES


- Predictive and physiological cell model; applicable for drug development and preclinical research (e.g., ventricular hypertrophy)
- Established assays for cardiac contractile force measurements, impedance measurements, calcium transients analysis and metabolic liability
- Ideal for bioengineering including disease and tissue modeling

# vCor.4U™

Human Induced Pluripotent Stem Cell-Derived Ventricular Cardiomyocytes




## AP STABILITY ENABLES vCOR.4U™ FOR PHARMACOLOGICAL STUDIES



**Figure 2: vCor.4U™ elicit stable action potentials**

(A) Representative trace showing evoked AP stability following 0.1% DMSO application (10 min). (B) Action potential duration stability of evoked AP over time (APD vs. time) during the control period and following application of 0.1% DMSO.

## vCOR.4U™ ARE SUITABLE TO ASSESS COMPOUNDS' PRO-ARRHYTHMIC RISK



**Figure 3: Pharmacological sensitivity of vCor.4U™ to cardiac channel modulators**

(A) Table summarizing the response of vCor.4U™ to key cardiac channel modulators. Carbachol was used to determine if the atrial current I<sub>KACH</sub> could be detected. No response was observed indicating a ventricular phenotype. (B) Application of dofetilide, an I<sub>Kr</sub> inhibitor with a high Torsade de Point (TdP) risk, elicited EADs (arrows) from spontaneous AP, suggesting that vCor.4U™ are a suitable model to assess the pro-arrhythmic risk of compounds.

## AXIOGENESIS OVERVIEW

### DIFFERENTIATED HUMAN CELLS

Axiogenesis is a leading expert in providing commercial-grade *in vitro* differentiated cell types derived from human induced pluripotent stem cells (iPSCs).

Core products include Cor.4U® cardiac myocytes and fibroblasts, as well as Peri.4U™, Dopa.4U™, CNS.4U™ and Astro.4U™ neural cells.

### VALIDATED ASSAYS & PROTOCOLS

Axiogenesis enables customer efficiency by providing ready to use cells along with validated protocols.

Assays for each cell type have been developed for advanced drug discovery, safety pharmacology, *in vitro* toxicology applications, and disease and tissue modeling.

Based on its in-house assay capabilities, Axiogenesis is able to provide expert scientific support in order to facilitate selection and quick implementation of validated assays and technologies.



iPSC-derived neurons



iPSC-derived cardiomyocytes

## PRODUCT SPECIFICATIONS

Cell type	iPSC-derived ventricular cardiomyocytes
Source	iPSC of 26 y/o Caucasian female
Species	Human
Purity	Approx. 76% ventricular; fibroblast-free
Assay window	Stable beating after 72h. Refer to our protocols for assay-specific recommendations



vCor.4U™



Metrion Ion Channel Screen Service



FOR MORE INFORMATION VISIT [WWW.AXIOGENESIS.COM](http://WWW.AXIOGENESIS.COM) OR CONTACT [INFO@AXIOGENESIS.COM](mailto:INFO@AXIOGENESIS.COM)

**Europe**  
Nattermannallee 1/S20  
50829 Cologne  
Germany  
+49 221 99 88 18-0

**North America**  
600 W Germantown Pike, STE 110  
Plymouth Meeting, PA 19462  
USA  
+1 844-511-6959

