



A Clinically Translatable hiPSC Cardiomyocyte Model to Predict QTc and QRS Cardiac Risk

Steve Jenkinson

24th September 2024

www.metrionbiosciences.com

hiPSC Cardiomyocytes in Defining CV Risk

Introduction

- hiPSC-CMs are now routinely used in industry to assess CV effects of novel compounds
 - Functional endpoints - Contractility
 - Electrophysiological endpoints – QTc/QRS risk assessment
- Importance noted by regulators resulting in change in ICH guidelines
 - Driven initially by CiPA initiative
 - Formalized by inclusion in revised ICH S7B Q&As
 - hiPSC-CM data can be used to support a Thorough QT (TQT) waiver application
- Key requirement of any hiPSC-CM assay is an understanding of its translation to the clinic

hiPSC Cardiomyocytes in Defining CV Risk

Introduction

- Presentation describes a high throughput hiPSC-CM model that can
 - Predict exposures of a novel compound that would be associated with a 10 ms change in clinical QTc interval
 - Define the probability of QRS prolongation risk
 - Assess general CV toxicity liability
 - Acutely (30 min) and chronic (24 h) endpoints examined
 - Assay performed in serum free conditions


hiPSC Cardiomyocytes in Defining CV Risk

Cell System

iCell² cardiomyocytes

(FUJIFILM Cellular Dynamics Inc.)

- Key ion channels involved in cardiac AP generation show similar expression to purified human ventricular cardiomyocytes




Kilfoil et al. (2021) *Eur. J. Pharmacol.* (PMID: 34678241)

hiPSC Cardiomyocytes in Defining CV Risk


The Technology

Action potential recordings from intrinsically paced hiPSC-CMs


- Volta Fast Optical Reader (Lumencor Inc.)
- Utilizes voltage sensitive membrane dye (BeRST)
- Simultaneous reads 96 well plate at 10,000 Hz
- Recordings equivalent in quality to patch clamp



Vehicle



AP Prolongation



Early After
Depolarizations




hiPSC Cardiomyocytes in Defining CV Risk

Assay Endpoints

Concentration-response curves for multiple endpoints


- 10-point concentration response curves
- Endpoints measured
 - Action potential duration (APD)
 - Rise time
 - Beat rate
- Acute (30 min) and chronic (24 h) measurements
- Protein/serum free assay




hiPSC Cardiomyocytes in Defining CV Risk

Example Data


hERG Inhibition



Ca_v1.2 Inhibition



hERG Trafficking Inhibition






Kilfoil et al. (2021) Eur. J. Pharmacol. (PMID: 34678241)

hiPSC Cardiomyocytes in Defining CV Risk

Rate Corrected APD₉₀

APD₉₀ rate correction

- Yamamoto correction (ycAPD₉₀) is optimal method for hiPSC-CMs
- ycAPD₉₀ is a surrogate of clinical QTc



Kilfoil et al. (2021) *Eur. J. Pharmacol.* (PMID: 34678241)

Predicting Clinical QTc Risk Using hiPSC-CMs



www.metribiosciences.com



hiPSC Cardiomyocytes in Defining CV Risk

Prediction of Clinical QTc

- Ability to determine QTc risk was assessed using a series of reference compounds
 - 23 QTc positive reference compounds with available clinical data
- The analysis compared
 - The concentration associated with hiPSC-CM APD₉₀ threshold value
 - 3x Vehicle Std Dev
 - Free clinical exposure producing 10 ms QTc change




Reference Compounds		
Azimilide*	Ibutilide*	Ribociclib
Bepidil*	Mesoridazine	Sotalol*
Cisapride*	Moxifloxacin	Terfenadine*
Citalopram	Odansetron*	Terodiline
Dofetilide*	Procainamide	Thioridazine
Droperidol	Quinidine*	Tolterodine
E4031	Quinine	Vandetanib*
Halofantrine	Ranolazine*	

*CiPA 28 compound

hiPSC Cardiomyocytes in Defining CV Risk

Prediction of Clinical QTc

Ability to predict the free clinical exposure associated with a 10 ms change in QTc for novel compounds



Input Threshold Value
Into Equation
($Y=0.92X - 0.35$)

Prediction of Free Clinical
Exposure Producing
a 10 ms QTc Change

- Predictivity confirmed using a 4-fold cross validation analysis
- hiPSC-CM data available for an additional 43 reference compounds (66 in total)
 - Including all CiPA-28 compounds

Assessing the Probability of Clinical QRS Risk Using hiPSC-CMs


www.mettrionbiosciences.com



hiPSC Cardiomyocytes in Defining CV Risk

Probability of Clinical QRS Liability

- Action Potential Rise Time is associated with $\text{Na}_v1.5$ channel activity
- $\text{Na}_v1.5$ block is associated with QRS prolongation in the clinic
- hiPSC-CM assay utilizes Rise Time to assess clinical QRS risk
- Available clinical QRS data less defined than for QTc
 - Clinical data used in analysis
 - Clinician defined effect on QRS (Binary - Yes/No)
 - Free clinical exposure associated with finding
 - Stem cell data used in analysis
 - Rise Time threshold concentration




hiPSC Cardiomyocytes in Defining CV Risk

Probability of Clinical QRS Liability

Performed ROC analysis to define hiPSC-CM/clinical data association

- 22 QRS positive / 12 QRS negative compounds
- Compared QRS effect (Yes/No) versus ratio of stem cell rise time threshold concentration versus free clinical exposure associated with finding



Incubation	AUROC (95% CI)	Sensitivity (95% CI)	Specificity (95% CI)	Cut Point	N Value
30 min	0.98 (0.95, 1.01)	0.92 (0.64, 1.00)	0.95 (0.76, 1.00)	33.0	34
24 h	0.98 (0.95, 1.02)	0.92 (0.64, 1.00)	0.95 (0.76, 1.00)	56.4	34

Kilfoil et al. (2021) *Eur. J. Pharmacol.* (PMID: 34678241)


hiPSC Cardiomyocytes in Defining CV Risk

Probability of Clinical QRS Liability

Defining QRS probability

- QRS probability at a clinical exposure equivalent to
 - hiPSC-CM Rise Time threshold concentration
 - 74 - 86 %
 - 33-fold (30 min) or 56-fold (24 h) lower than Rise Time threshold concentration
 - 1.2 - 2.4%
 - Analysis assumes Nav1.5 prevalence of
 - 12.7 - 23.5%

Impact of hiPSC-CM Data On Defining QRS Risk



Case Study

www.mettrionbiosciences.com



hiPSC Cardiomyocytes in Defining CV Risk

Case Study

- Compound ion channel profile
 - hERG $IC_{50} = \sim 20 \mu\text{M}$
 - $Na_v1.5$ $IC_{50} = \sim 115 \mu\text{M}$
 - $Ca_v1.2$ $IC_{50} = \sim 280 \mu\text{M}$
- QTc and QRS prolongation observed in the clinic
 - Free clinical exposure $\sim 1 \mu\text{M}$
- QRS liability was unexpected based on $Na_v1.5$ patch clamp data
 - hiPSC predicted 10 ms change in QTc plus QRS risk around $1 \mu\text{M}$ free exposure
- Further analysis suggests that $Na_v1.5$ block had components of state, use/rate dependency
- **Holistic hiPSC-CM model was able to predict effects missed in routine ion channel profiling**

hiPSC Cardiomyocytes in Defining CV Risk

Summary

- hiPSC-CM model is valuable in assessing CV risk
 - Ability to predict exposure associated with a 10 ms QTc change in the clinic
 - Assess probability of clinical QRS risk and associated exposure
 - Also valuable in assessing long term general toxicity
- Provide insight into potential mechanism of CV effects (hERG, Na_v1.5, Ca_v1.2)
- Assessment of acute (30 min) or chronic (24 h) effects of novel compounds
- Assay is valuable in assessing CV liabilities of novel compounds in a holistic integrated system

Questions?

Visit us at Booth 203

Contact us:

Steve Jenkinson

VP Drug Discovery and Safety

steve.jenkinson@metrionbiosciences.com

www.metrionbiosciences.com




Backup Slides

www.mettrionbiosciences.com



Potential Profiling Scheme

Number of
Compounds
Profiled



hERG Profiling



Na_v1.5/Ca_v1.2 Profiling



hiPSC-CM Assay
(Ion channel Plus Overt
Acute or Chronic Tox Derisking)

Effects Observed

CiPA Ion Channel Panel
(To rule out ion channel effects)



In Vivo CV Study



GLP hERG
(IND-Enabling Study)