

## Seal enhancers on the Qube 384: An alternative to F<sup>-</sup>

Catherine Hodgson<sup>1</sup>, Edward Stevens<sup>1</sup>, Anissa Bara<sup>2</sup>, Beatrice Badone<sup>2</sup>, Gary Clark<sup>1</sup>

<sup>1</sup>Metrion Biosciences Ltd, Granta Centre, Granta Park, Cambridge, CB21 6AL, UK

<sup>2</sup>Sophion Bioscience A/S, Industriparken 39, 2750 Ballerup, Denmark

### Introduction

Gigaohm seals, or 'gigaseals', are imperative to patch clamp electrophysiology to enable good electrical access to the cell and high-quality recordings. These seals form through chemical bonds and electrostatic forces between the cell membrane and the glass pipette in manual patch clamp, or in the case of planar patch clamp, between the cell membrane and chip substrate<sup>1</sup>. Planar patch clamp often requires the use of 'seal enhancers' to increase the resistances of these seals, with CaF<sub>2</sub> being the most extensively used. It is hypothesised that high concentrations of extracellular Ca<sup>2+</sup> and intracellular F<sup>-</sup> give rise to CaF<sub>2</sub> precipitate at the solution interface, fostering seal formation<sup>2</sup>.

CaF<sub>2</sub> as a seal enhancer, however, has limitations. F<sup>-</sup> is known to stimulate G-protein modulation of ion channels, altering channel properties<sup>3-6</sup>. Furthermore, use of F<sup>-</sup> is not optimal when recording from Ca<sup>2+</sup>-activated ion channels due to resultant unknown concentrations of free intracellular Ca<sup>2+</sup>.

In an effort to overcome these limitations, Metrion and Sophion collaborated to determine whether other insoluble salts can act as seal enhancers, or whether this property is unique to CaF<sub>2</sub>.


### Results

The solubility product constants ( $K_{sp}$ ) of different Ca<sup>2+</sup>, Ba<sup>2+</sup> and Sr<sup>2+</sup> salts are compared in Table 1a. BaSO<sub>4</sub> and SrCO<sub>3</sub> have low  $K_{sp}$  values (with corresponding low solubility), similar to CaF<sub>2</sub>. The PO<sub>4</sub><sup>3-</sup> salts, however, have even lower  $K_{sp}$  values


and are highly insoluble. The ability of the different salt pairs to promote gigaseal formation was analysed using a CHO-hNav1.5 cell line and single hole QChips. The cells were resuspended in extracellular solutions containing 10 mM Ca<sup>2+</sup>, Ba<sup>2+</sup> or Sr<sup>2+</sup>, whilst a selection of intracellular solutions was used containing high concentrations of the different anions. The concentration of each anion was adjusted to ensure consistent osmolarity across the intracellular solutions (~290 mOsm). Only BaSO<sub>4</sub> successfully facilitated the formation of gigaseals comparable to CaF<sub>2</sub> (Table 1b). Interestingly, these data suggest there is no correlation between salt  $K_{sp}$  values and their ability to foster gigaseal formation, particularly as neither SrCO<sub>3</sub> nor the PO<sub>4</sub><sup>3-</sup> salts significantly increased seal resistances. Thus, the relationship between salt crystal formation and gigaseals remains unclear. Still, as SO<sub>4</sub><sup>2-</sup> is not known to alter ion channel function or cellular signalling pathways, BaSO<sub>4</sub> represents a candidate alternative seal enhancer to CaF<sub>2</sub> with potential advantages.

Experiments were subsequently conducted to determine the minimal concentration of extracellular Ba<sup>2+</sup> required for the formation of gigaseals. To this end, CHO-hNav1.5 cells were resuspended in extracellular solutions containing different concentrations of Ba<sup>2+</sup>, and seals formed in the presence of different levels of intracellular SO<sub>4</sub><sup>2-</sup>. High resistance gigaohm seals only formed with ≥ 3 mM extracellular Ba<sup>2+</sup> (Table 2). It is known, however, that Ba<sup>2+</sup> ions block K<sup>+</sup> channels in the μM-mM range<sup>7-9</sup>. Thus, whether this concentration is low enough to spare K<sup>+</sup> channels from block is yet to be determined.

**Table 1. Correlation between salt pair solubility product constants ( $K_{sp}$ ) and gigaseal formation on Qube 384.** No correlation was found between the solubility product constants ( $K_{sp}$ ) of Ca<sup>2+</sup>, Ba<sup>2+</sup> and Sr<sup>2+</sup> salts (A) and their ability to foster gigaseal formation (B). Only BaSO<sub>4</sub> successfully facilitated the formation of gigaseals comparable to CaF<sub>2</sub>. Despite the PO<sub>4</sub><sup>3-</sup> salts having very low  $K_{sp}$  values and SrCO<sub>3</sub> having a similar  $K_{sp}$  value to CaF<sub>2</sub> and BaSO<sub>4</sub>, these salts failed to produce gigaohm seals. Moderate seal resistances with PO<sub>4</sub><sup>3-</sup> salts were transient and unstable. BaF<sub>2</sub> and SrF<sub>2</sub> were not studied as alternative seal enhancers to CaF<sub>2</sub> as they do not eliminate the use of F<sup>-</sup>. Median resistances calculated from 24 cells per salt pair.



**Table 2. Resistances of seals formed using various concentrations of extracellular Ba<sup>2+</sup> and intracellular SO<sub>4</sub><sup>2-</sup>.** Gigaseals only formed with ≥3 mM Ba<sup>2+</sup> and in the presence of SO<sub>4</sub><sup>2-</sup>. Median resistances calculated from 24 or 48 cells per condition.




## Effects of CaF<sub>2</sub> and BaSO<sub>4</sub> on hNa<sub>v</sub>1.5 biophysical properties

F is known to trigger G-protein modulation of ion channels and has been reported to affect the voltage-dependent kinetics of Na<sub>v</sub> channels<sup>4-6, 10-13</sup>. Single hole experiments were conducted to assess the effects of CaF<sub>2</sub> and BaSO<sub>4</sub> as seal enhancers on hNa<sub>v</sub>1.5 biophysical properties.

CHO-hNa<sub>v</sub>1.5 cells were resuspended in extracellular solutions containing 3, 5 or 10 mM Ca<sup>2+</sup> or Ba<sup>2+</sup>, and seals formed in the presence of intracellular solutions containing various concentrations of F<sup>-</sup> or SO<sub>4</sub><sup>2-</sup>. A series of 3-s pulses between -100 mV and +50 mV in +10 mV intervals from a holding potential of -120 mV, followed by a 20-ms pulse at 0 mV was used to determine the voltage dependence of activation and inactivation of hNa<sub>v</sub>1.5. There was a significant depolarising shift in hNa<sub>v</sub>1.5 V<sub>0.5</sub>

inactivation with increasing concentrations of intracellular F<sup>-</sup> (Figure 1a). In contrast, increasing SO<sub>4</sub><sup>2-</sup> concentration had no significant effect on V<sub>0.5</sub> inactivation (Figure 1b). This suggests that, unlike the G-protein effects of F<sup>-</sup> on hNa<sub>v</sub>1.5, SO<sub>4</sub><sup>2-</sup> has no effect on intracellular signalling and does not impact channel biophysics. Of note, there was also a depolarising shift in hNa<sub>v</sub>1.5 V<sub>0.5</sub> inactivation with higher concentrations of extracellular Ca<sup>2+</sup> and Ba<sup>2+</sup> most likely due to surface potential screening effects (independent of G-protein modulation)<sup>14-15</sup>. In contrast to the voltage dependence of inactivation, there was no significant difference in hNa<sub>v</sub>1.5 current amplitude or V<sub>0.5</sub> activation between the two seal enhancers at the different concentrations (data not shown).




**Figure 1. Effects of CaF<sub>2</sub> and BaSO<sub>4</sub> on hNa<sub>v</sub>1.5 channel biophysics.** hNa<sub>v</sub>1.5 V<sub>0.5</sub> inactivation with different cation and anion concentrations (mean ± S.D., N ≥ 11). Increasing concentrations of intracellular F<sup>-</sup> caused a depolarising shift in V<sub>0.5</sub> inactivation (A). In contrast, increasing concentrations of SO<sub>4</sub><sup>2-</sup> had no effect on hNa<sub>v</sub>1.5 V<sub>0.5</sub> inactivation (B). One-way ANOVAs conducted within each cation group followed by Tukey's Honestly Significant Difference post-hoc tests: \*\* = *p* < .01; \*\*\* = *p* < .001; \*\*\*\* = *p* < .0001. CaF<sub>2</sub> and BaSO<sub>4</sub> data taken from two separate experiments.

## hNav1.5 pharmacology: CaF<sub>2</sub> versus BaSO<sub>4</sub>

Next, multihole Qube 384 experiments were conducted to assess differences in hNav1.5

pharmacology between CaF<sub>2</sub> and BaSO<sub>4</sub>. Potency IC<sub>50</sub> values for a range of Na<sup>+</sup> channel inhibitors were comparable against hNav1.5 between the different seal enhancers (Figure 2).




**Figure 2. CaF<sub>2</sub> versus BaSO<sub>4</sub> – hNav<sub>v</sub>1.5 pharmacology.** **A)** Representative sweep plots (left) and current-time (I-t) plots (right) for hNav<sub>v</sub>1.5 inhibition by amitriptyline. There was no difference in cumulative inhibition of hNav<sub>v</sub>1.5 by increasing concentrations of amitriptyline between CaF<sub>2</sub> and BaSO<sub>4</sub>. **B)** Screening of a range of inhibitory compounds showed no difference in hNav<sub>v</sub>1.5 pharmacology between CaF<sub>2</sub> and BaSO<sub>4</sub>. Concentration-response curves for amitriptyline against hNav<sub>v</sub>1.5 using CaF<sub>2</sub> or BaSO<sub>4</sub> as the seal enhancer (mean ± S.D., N = 12 wells per concentration for CaF<sub>2</sub>, N = 8 wells per concentration for BaSO<sub>4</sub>). CaF<sub>2</sub> and BaSO<sub>4</sub> data taken from two separate experiments.


### hCa<sub>v</sub>1.2 pharmacology: CaF<sub>2</sub> versus BaSO<sub>4</sub>

In a similar manner, the effects of CaF<sub>2</sub> versus BaSO<sub>4</sub> on hCa<sub>v</sub>1.2 pharmacology were assessed using a HEK293-hCa<sub>v</sub>1.2 cell line and multihole QChips. Using BaSO<sub>4</sub> as the seal enhancer changes the kinetics of the channel (Figure 3). Ca<sup>2+</sup> influx with CaF<sub>2</sub> activates calmodulin, which binds to the intracellular regions of the channel, facilitating Ca<sup>2+</sup>-dependent inactivation<sup>16-17</sup>. Ba<sup>2+</sup> is routinely used as a surrogate charge carrier in the study

of Ca<sup>2+</sup> channels<sup>18</sup>. As expected, Ca<sup>2+</sup>-dependent inactivation of hCa<sub>v</sub>1.2 was eliminated when Ba<sup>2+</sup> was used as a surrogate carrier ion with BaSO<sub>4</sub> as the seal enhancer, also giving rise to larger hCa<sub>v</sub>1.2 current amplitudes (mean ± S.D. – CaF<sub>2</sub>: 0.79 ± 0.58 nA (N = 11) versus BaSO<sub>4</sub>: 1.29 ± 0.89 nA (N = 56); Welch's *t*-test: *t*<sub>(20)</sub> = 2.36; *p* < .05). Despite this, hCa<sub>v</sub>1.2 pharmacology assessed using two inhibitory compounds, nifedipine and verapamil, was unaffected (Figure 4).



**Figure 3. CaF<sub>2</sub> versus BaSO<sub>4</sub> – hCa<sub>v</sub>1.2 kinetics.** hCa<sub>v</sub>1.2 exhibits Ca<sup>2+</sup>-dependent inactivation when CaF<sub>2</sub> is used as the seal enhancer (A). BaSO<sub>4</sub> as the seal enhancer (using Ba<sup>2+</sup> as a surrogate carrier ion) confers loss of the Ca<sup>2+</sup>-dependent inactivation of hCa<sub>v</sub>1.2 observed with CaF<sub>2</sub> (B). Example sweep plots derived from Sophion Analyzer v9.0.42.




**Figure 4. CaF<sub>2</sub> versus BaSO<sub>4</sub> – hCa<sub>v</sub>1.2 pharmacology.** Mean ± S.D. concentration-response curves for two common inhibitors against hCa<sub>v</sub>1.2, nifedipine (A) and verapamil (B) (CaF<sub>2</sub>: N = 2-5 wells per concentration; BaSO<sub>4</sub>: N = 6-12 wells per concentration). Compound potencies (IC<sub>50</sub> values) did not differ between CaF<sub>2</sub> and BaSO<sub>4</sub>.

## Seal stability following removal of seal enhancers

Exchange of solutions following gigaseal formation would provide flexibility to assay different ion channels using optimal recording solutions. However, stability of gigaseals formed in the presence of BaSO<sub>4</sub> (3, 5 or 10 mM

extracellular Ba<sup>2+</sup> and 50 or 70 mM intracellular SO<sub>4</sub><sup>2-</sup>) were not maintained upon exchange for corresponding concentrations of extracellular Ca<sup>2+</sup> (Figure 5), whilst those formed using CaF<sub>2</sub> as the seal enhancer (3, 5 or 10 mM extracellular Ca<sup>2+</sup> and 120 mM intracellular F<sup>-</sup>) were lost or greatly reduced following exchange for intracellular Cl<sup>-</sup> across a range of concentrations (data not shown).



**Figure 5. Ba<sup>2+</sup>-Ca<sup>2+</sup> solution exchange.** Example  $R_{\text{membrane}}$ -time plot showing loss of a gigaohm seal formed using 5 mM extracellular Ba<sup>2+</sup> and 70 mM intracellular SO<sub>4</sub><sup>2-</sup> upon exchange for 5 mM extracellular Ca<sup>2+</sup>.

## Conclusion

In conclusion, although no correlation was found between salt  $K_{\text{sp}}$  values and gigaseal formation, BaSO<sub>4</sub> was identified as an equivalent seal enhancer to CaF<sub>2</sub> in planar patch clamp electrophysiology. BaSO<sub>4</sub> and CaF<sub>2</sub> were characterised across two different ion channels, hNav1.5 and hCa<sub>v</sub>1.2. Whilst increasing concentrations of F<sup>-</sup> caused depolarising shifts in the voltage dependence of inactivation of hNav1.5, SO<sub>4</sub><sup>2-</sup> had no effects on hNav1.5 biophysical properties. Additionally, there was no difference in the pharmacological effects of inhibitory compounds against hNav1.5 or hCa<sub>v</sub>1.2 between CaF<sub>2</sub> and BaSO<sub>4</sub> seal enhancers. BaSO<sub>4</sub> could be used as the seal enhancer when recording from non-K<sup>+</sup>-conducting Ca<sup>2+</sup>-activated channels, such as the Ca<sup>2+</sup>-activated Cl<sup>-</sup> channel TMEM16A, allowing more accurate estimation of the free intracellular Ca<sup>2+</sup> concentration.

## Methods

Experiments were conducted using a Sophion Bioscience Qube 384 with QChip 384 (single hole) and QChip 384X (multihole) consumables. Temperature was maintained at 22 °C using the Qube temperature control module.

Analysis was conducted using Sophion Analyzer v9.0.42 and GraphPad Prism v10.2.2.

CHO-hNav1.5 and HEK293-hCa<sub>v</sub>1.2 cell lines were provided by Metrion Biosciences.

All compounds were tested at: 0.001, 0.01, 0.1, 1, 10 and 100 μM.

## Acknowledgements

Thank you to Andrew James for ideas and edits, and to the Metrion Biosciences cell culture team, Aqila March, Megan Pritchard, Jack O'Donnell and Joana Rocha.

## References

- <sup>1</sup> Suchyna, T. M., Markin, V. S., & Sachs, F. (2009). Biophysics and structure of the patch and the giga-seal. *Biophysical Journal*, 97(3), 738-747.
- <sup>2</sup> Løjknær, L. D., Lindqvist, A., Sauter, D. R., & Rosholm, K. R. (2023). Seal enhancer. US 11692994 B2.
- <sup>3</sup> Vargas, G., Yeh, T. Y. J., Blumenthal, D. K., & Lucero, M. T. (1999). Common components of patch-clamp internal recording solutions can significantly affect protein kinase A activity. *Brain Research*, 828(1-2), 169-173.
- <sup>4</sup> Rapedius, M., Obergrussberger, A., Humphries, E. S. A., Scholz, S., Rinke-Weiss, I., Goetze, T. A., Brinkwirth, N., Rotordam, M. G., Strassmaier, T., Randolph, A., Friis, S., Liutkute, A., Seibertz, F., Voigt, N., & Fertig, N. (2022). There is no F in APC: using physiological fluoride-free solutions for high throughput automated patch clamp experiments. *Frontiers in Molecular Neuroscience*, 15, 982316.
- <sup>5</sup> Antonny, A., Sukumar, M., Bigay, J., Chabre, M., & Higashijima, T. (1993). The mechanism of aluminum-independent G-protein activation by fluoride and magnesium. 31P NMR spectroscopy and fluorescence kinetic studies. *Journal of Biological Chemistry*, 268(4):2393-2402.
- <sup>6</sup> Vincent, S., Brouns, M., Hart, M. J., & Settleman, J. (1998). Evidence for distinct mechanisms of transition state stabilization of GTPases by fluoride. *Proceedings of the National Academy of Sciences of the United States of America*, 95(5):2210-2215.
- <sup>7</sup> Armstrong, C. M., Swenson, R. P., & Taylor, S. R. (1982). Block of squid axon K channels by internally and externally applied barium ions. *Journal of General Physiology*, 80(5), 663-682.
- <sup>8</sup> Tokimasa, T. & Nishimura, T. (1997). Actions of barium on rapidly inactivating potassium current in bullfrog sympathetic neurons. *Neuroscience Letters*, 236(1), 37-40.
- <sup>9</sup> Kotani, S., Hirasawa, T., Suzuki, T., Sato, K., Sakakibara, M., & Tokimasa, T. (2000). Mechanisms underlying the M-current block by barium in bullfrog sympathetic neurons. *Neuroscience Letters*, 285(1), 1-4.
- <sup>10</sup> Meadows, L. S., Chen, Y. H., Powell, A. J., Clare, J. J., & Ragsdale, D. S. (2002). Functional modulation of human brain Nav1.3 sodium channels, expressed in mammalian cells, by auxiliary beta 1, beta 2 and beta 3 subunits. *Neuroscience*, 114(3), 745-753.
- <sup>11</sup> Rugiero, F., Mistry, M., Sage, D., Black, J. A., Waxman, S. G., Crest, M., Clerc, N., Delmas, P., & Gola, M. (2003). Selective expression of a persistent tetrodotoxin-resistant Na<sup>+</sup> current and Nav1.9 subunit in myenteric sensory neurons. *The Journal of Neuroscience*, 23(7), 2715-2725.
- <sup>12</sup> Coste B., Osorio, N., Padilla, F., Crest, M., & Delmas, P. (2004). Gating and modulation of presumptive Nav1.9 channels in enteric and spinal sensory neurons. *Molecular and Cellular Neuroscience*, 26(1), 123-134.
- <sup>13</sup> Jarecki B. W., Jackson J. O., II, Cummins T. R., Sheets P. L. (2008). Paroxysmal extreme pain disorder mutations within the D3/S4-S5 linker of Nav1.7 cause moderate destabilization of fast inactivation. *The Journal of Physiology*, 586(17), 4137-4153.
- <sup>14</sup> Hille, B., Woodhull, M., & Shapiro, B. I. (1975). Negative surface charge near sodium channels of nerve: divalent ions, monovalent ions, and pH. *Philosophical Transactions of the Royal Society B: Biological Sciences*. 270(908), 301-318.
- <sup>15</sup> Hahn, R. & Campbell, D. T. (1983). Simple shifts in the voltage dependence of sodium channel gating caused by divalent cations. *Journal of General Physiology*. 82(6), 785-805.
- <sup>16</sup> Kameyama M., Minobe, E., Shao, D., Xu J., Gao, Q., & Hao, L. (2023). Regulation of cardiac Cav1.2 channels by calmodulin. *International Journal of Molecular Sciences*, 24(7), 6409.
- <sup>17</sup> Sang, L., Vieira, D. C. O., Yue, D. T., Manu, B. J., & Dick, I. E. (2021). The molecular basis of the inhibition of Cav1 calcium-dependent inactivation by the distal carboxy tail. *Journal of Biological Chemistry*, 296, 100502.
- <sup>18</sup> Ren, M., Randolph, A. L., Alvarez-Baron, C., Guo, D., Tran, P. N., Thiebaud, N., Sheng, J., Zhao, J., & Wu, W. (2022). Experimental factors that impact Cav1.2 channel pharmacology – effects of recording temperature, charge carrier, and quantification of drug effects on the step and ramp currents elicited by the “step-step-ramp” voltage protocol. *PLoS One*, 17(11), e0276995.